

Dermo Jet

Intradermal Needle-less Injector

Injection without needle provides obvious advantages:-

- Relatively painless for the patient.
- Cost saving in terms of needles, syringes and medication.
- Intradermal injections of 0.1ml per injection are achieved by a high velocity liquid jet which penetrates the surface membrane to a depth of 2-5mm and produces a weal of up to 6mm diameter without contact with the injection site.
- The 4ml reservoir is easily filled to permit up to 40 injections to be given in rapid succession.

Cat No. JI-130


Specifications:

| | |
|----------------------------|---------------------------------------|
| <i>Dosage:</i> | <i>0.1 ml per injection</i> |
| <i>Reservoir Capacity:</i> | <i>4 ml</i> |
| <i>Loading:</i> | <i>By hand operated lever</i> |
| <i>Weight:</i> | <i>300 grams</i> |
| <i>Construction:</i> | <i>Surgical grade stainless steel</i> |

Fields of Application:

- | | |
|---------------------------------|---------------|
| ■ Dermatology & Plastic Surgery | ■ Anaesthesia |
| ■ Accident & Emergency | ■ Vaccination |
| ■ Pain Clinics | ■ Dentistry |
| ■ Orthopaedic Surgery | ■ Podiatry |

Cyto-Solutions Ltd



Tel: 01233 665522

Fax: 01233 665523

Email: info@cyto-solutions.co.uk